

# Ballast Water test kit

## Complete Testing - No Extra Costs!

A complete Portable Ballast Kit, according to **IMO regulations**, the Unimarine Ballast Water Test Kit includes the *Vibrio cholera test* and the *Enterococci tests* ("Enterolert").

The **Ballast Water Test Kit** consist of two cases:

### 1 Ballast Water Test Kit #1 Manual testing process

- Chloride Titration Test
- Aquavial E. Coli & Coliform
- Vibrio Cholerae O1 Rapid Test
- Enterolert Test Kit



### 2 Ballast Test Kit #2 Handheld Indicative Monitoring test kit

Light, highly durable Fluorometer.  
Checks the Algae (Chlorophyll)  
& Plankton in water.



Connect with us!

[www.unimarinegroup.com](http://www.unimarinegroup.com)



@unimarinegroup

+357 25331054

+357 99833991 (24/7)

[info@unimarine-services.com](mailto:info@unimarine-services.com)

## Contents & Specifications:

### 1 Ballast Water Test Kit #1 Manual testing process

#### Contents:

- **Chloride Titration Test**  
100 ml beaker  
Crushing and stirring rod  
Cleaning brush  
Pack with 100 Chloride TC tablets

**Measurement Range:** 0 – 5000 mg/l Cl<sup>-</sup>  
0 – 8000 mg/l NaCl

- **Aquavial E. Coli & Coliform**  
Plastic casing  
User manual  
4 x 5 ml plastic vials  
with culture medium

Detects E. Coli / Coliform concentrations as low as 1 CFU/ml (one bacteria per ml) of E. Coli / Coliform bacteria within **24 hours** when incubated at: 35 – 40 °C (95 – 104 °F), or **72 hours** at: 20 – 25 °C (68 – 77 °F).

- **Vibrio Cholerae O1 Rapid Test**  
Foil wrapped test cassettes with dropper desiccant (x10)  
Faecal collection device containing buffer (x10)

(Positive / Negative)

- **Entolert Test Kit**  
Enterolert Enterococci tests\*  
Thio-Bags (x25)

Detects enterococci, such as E. faecium and E. faecalis, in fresh & marine water.

\*pack of 10

### 2 Ballast Test Kit #2 Handheld Indicative Monitoring test kit

#### Contents:

- Ballast Check-2 Flurometer
- Filter Washing Kit
- Filter Capsules (x3)
- 50ml Syringe
- Glass Cuvette
- Box of AAA4 batteries
- Flash Drive (USB)
- Data Download Cable
- Wrist Strap
- Quick Start Guide




Connect with us!

[www.unimarinegroup.com](http://www.unimarinegroup.com)



@unimarinegroup   

 +357 25331054

 +357 99833991 (24/7)

 [info@unimarine-services.com](mailto:info@unimarine-services.com)